

StandingVest


SystemRoMedic™

Trial Fitting Protocol

Patient's name: _____

Social insurance number: _____

StandingVest is specially designed for lifting with a standing shell. For the caregiver, it is easy to apply and it gives the patient a feeling of security. It is designed to provide upper-body support and is held in place by leg harness at the groin. StandingVest is designed to activate the patient safely and comfortably.



Functional inspection

Visual and mechanical inspection

Check the condition and function of the sling regularly. Always inspect the product after laundering. Check to ensure that seams and material are free from damage. Check to ensure that fabric is not worn or faded. Apply load to the device and check to ensure that clasps, handles, etc. withstand heavy load. If there are signs of wear, the product must be discarded.

Which size is used?

X-Small Small Medium Large X-Large

Which material is used?

Polyester Other

Accessories

SlingBarWrap AdjustmentLoop StrapsPadding

What is the width of the slingbar?

Instructions for use

Can the sling/vest be used in showering situations?

Yes No

AdjustmentLoop: _____

Have staff that use the assistive device:

Read the instructions?

Yes No

Received practical instruction in patient transfers?

Yes No

Instructions to personnel:

e.g.; Placement of head/ arms/hands? Placement of legs/feet? Use bed functions.

Number of caregivers: _____

Instructions to patient:

e.g.; In what way can the patient participate actively in the transfer situation? Which instructions are used?

Other:

Trial-fitted by:

Date: _____

Phone: _____

Phone: _____

Handicare AB
Veddestav. 15, Box 640
SE-175 27 Järfälla
SWEDEN

Tel: +46 (0)8-557 62 200
Fax: +46 (0)8-557 62 299
E-mail: info@handicare.se

Internet: www.handicare.com

In the events of changes or questions, contact the trial-fitter.

Manuals and instructions for functional testing are
available on Handicare's website www.handicare.com

